

Laboratory 9 – Muscle Physiology

Purpose:




The purpose of laboratory 9- Muscle Physiology, was to perform several experimental procedures in order to investigate the different contraction characteristics of skeletal, cardiac, as well as smooth muscle. In this lab there was a demonstration of the electromyograph that demonstrated the concepts of agonist, antagonist, and synergist muscles. The lab also demonstrated the effects of oxygen supply on skeletal muscle activity.

Procedure:

9-D: Demonstration of the electromyograph (EMG)

- Connect IWX/214 unit and open Labscribe3 program to generate electromyograph
- The subject must remove all jewelry and clean their skin that will be used. After the area dries place the electrode to the six locations.
- Place electrodes from proximal to distal on the forearm in this order: +2, -2 on posterior and +1, -1 and ground on the anterior. Snap the lead wires onto the electrodes in this order: red “+1” lead is attached to the proximal electrode on the anterior surface. Black “-1” lead is attached to the distal electrode on the anterior forearm. Green “C” lead (the ground) is attached to the remaining electrode on the anterior surface. White “+2” lead is attached to the proximal electrode on the posterior forearm. Brown “-2” lead is attached to the distal electrode on the posterior surface.
- Record an EMG of the muscles of the forearm that illustrates agonistic and antagonistic muscle activity for each exercise.
 - a. Flex the wrist with palm open and hold for four seconds. Return wrist to neutral position. Extend the wrist and hold for four seconds several times.
 - b. Forcefully flex the wrist with hand closed into a fist for four seconds. Return to neutral position. Extend the wrist maintaining a fist for four seconds several times.
 - c. Flex wrist against resistance applied by another student for 10 seconds.
 - d. Place the hand in mid-supination and make a fist. Move hand upwards against resistance applied by another student for 10 seconds

Results:



Discussion:

In this lab, different phases of typical skeletal muscle twitch were identified by performing these different exercises. The exercise that I felt demonstrated the most was pressing/flexing my wrist against another student because it showed more activity happening and it took more strength/power to try and go against the resistance.

Conclusion:

In conclusion, I believe this lab achieved its goal of demonstrating the different aspects of how muscle twitching works as well as being able to distinguish the EMG of a partially contracted muscle from a fully contracted muscle. I definitely learned a lot from this lab and what the EMG represents and how it demonstrates different aspects.